

# The V-Gen™ Dew Point / RH Generator

*Designed, built, and supported by InstruQuest Inc.*

## Model 1: Manual operation and interface for OEM implementation of automatic control

The V-Gen is a novel and versatile water vapor generation and calibration system. Its low cost, high quality and performance, and simplicity of operation make it ideally suited for a wide variety of applications requiring an accurate and stable source of water vapor. Its modular design, easy upgrade path, and versatility of control and software of higher models allow easy integration into a larger analytical setup or interfacing with other devices.

The V-Gen uses the fundamental, Two-Temperature principle to generate stable dew points suitable for calibration from 0 to 80 °C. A novel design substantially reduces stabilization time when changing the dew point temperature. Additional, Divided-Flow method is employed to generate low values of relative humidity (RH) that are not available by the Two-Temperature method alone. Thus, RH values from 0 to 100% are readily available up to 80 °C.

A precision pressure regulator isolates the selectable flow rate from fluctuations in the supply line. An automatic water supply system eliminates frequent user intervention inherent to other generators. The generated water vapor stream is delivered to the location of choice via flexible heated transfer line. To avoid condensation problems, the temperature of the transfer line is automatically maintained at a higher temperature than the current dew point value

## Applications

Stable Dew Point / RH source suitable for:

- Gravimetric sorption systems
- Permeation measurement instruments
- Relative humidity sensors calibration
- Sample conditioning in specialized instrumentation
- Study of materials hydration




## Operation

A dry or semi-dry gas source is connected at the instrument back via standard 1/8" compression fitting. Typically, the gas source can be air generated from an oil-less compressor or an inert gas can be supplied from a gas cylinder. The dryness of the gas will affect the lowest achievable RH values.

The high quality gas pressure regulator allows the selection of desired pressure and observing it on the gauge. Flow rate can be selected by adjusting the metering valve of the rotameter. Depending on application, rotameters with different ranges can be used and replaced by user.

The rotameter output is connected to the 3-valve manifold. In the two-temperature mode only the "WET" valve is open. If dry purge is desired at the selected flow rate, only the "DRY" valve is open. For divided flow operation the "WET" valve is open and the proportional valve, PV, is used. By increasing voltage to the PV driver, the



amount of flow bypassing the saturator increases. Thus, while maintaining the total flow constant, different ratios of dry gas flow to the saturated with water vapor gas flow can be obtained. An external RH sensor at the target needs to be used to monitor the output stream when the divided-flow method is employed.

The gas passing through the “WET” valve enters the saturator. Combined techniques of fine gas dispersion and passing it through media of large surface area are used to assure the highest degree of saturation for large range of temperatures and flow rates. Since evaporation of water causes small lowering of temperature inside the saturator and to account for any inefficiency, a second chamber called condenser is maintained at lower temperature to assure the 100% saturation (dew point). A RTD sensor is used for enhanced accuracy and resolution of the condenser temperature monitoring. The Watlow temperature controllers maintain the saturator and condenser temperatures selected by user.

The saturated stream leaving condenser enters the mixing manifold where it can be combined with dry gas. From the manifold the output stream is delivered via the heated transfer line to the location of choice. The temperature of the transfer line is automatically maintained at saturator temperature by a special controller board. Since the saturator temperature is normally maintained a few degrees above the condenser temperature, such setup eliminates any condensation problems.

The water used for gas saturation should be de-ionized and of high purity to avoid clogging of gas sparger and to minimize maintenance. The autonomous water delivery system consists of peristaltic pump and its controller, water level sensor, visual level indicator, pump enable switch (ON/OFF), and pump direction switch (Forward/Reverse) located at the back. When water level drops below certain level, the pump is activated and it supplies enough water to restore the required level. This amount of water is preheated to the saturator temperature to avoid any temperature drops in the saturator chamber. If automatic refill action is temporarily not desired, the pump enable switch can be turned off. In a case the water has to be removed from the system, e.g. for shipping purposes, the pump direction switch has to be switched to the REV position until the water is emptied.

No software is required to operate this model, but to facilitate calculations of the required dew point temperatures to obtain the relative humidity at the target location, a Windows® based software utility – “RH Calculator” is provided. Also, variety of other useful RH calculations can be carried out too.

## Control


The modular and “open-architecture” approach to control system allows easy interfacing to a larger analytical setup and addressing specific RH needs without unnecessary cost and complexity. In addition to manual mode of operation present in every model, the implementation of automatic control defines different models. For OEM users who prefer having direct control over hardware, all relevant digital and analog signals lines are connected to the DB15 socket at the instrument back. Model 1 also has additional RTD circuit for precise determination of dew point temperatures. Other models feature microcontroller board and ancillary electronics, serial communication port(s) and communication via a special command language, and PC software for experiment design and control.

## Calibration

One RTD used for precise dew point determination is calibrated using NIST-traceable precision glass thermometers. The second RTD used in control of condenser temperature is calibrated in situ using an independent RTD that has been calibrated using NIST-traceable precision glass thermometers. The calibration certificate for both RTD's is provided.


## System Dynamics

One of the objectives of V-Gen design was to improve the two-temperature system dynamics while assuring accuracy and stability.



The figure shows the heating and cooling rate of the condenser chamber. The whole cycle of heating the chamber and cooling it back takes less than an hour for the 80°C temperature span.

V-Gen temperature steps equilibration



This figure shows equilibration dynamics for several large temperature steps in the increasing and decreasing direction versus time. The room temperature was varied often throughout the experiment by about 2 °C.

## Specifications

### Dew Point Temperature:

Range: 0 to 90 °C  
 Accuracy: ± 0.2°C  
 Stability: 0.1C  
 Settability: 0.1°C

**RH range:** 0 to 100 %

**Maximum temperature of saturator:** 95 °C

**Maximum temperature of condenser:** 90 °C

**Thermal protection:** Thermal cut-offs for the saturator and condenser heaters. Temperature limits are also set in the temperature controllers.

**Flow Rate:** Determined by the rotameter range,  
 Available ranges: 50, 100, 200, 500, 1000 cc/min.

**Gas Type:** Inert gas, typically air or N<sub>2</sub>

**Gas Inlet Port:** 1/8" (Swagelok® type bulkhead)

### Gas Inlet Pressure:

Maximum: 100 psi (7 bar)  
 Minimum: 10 psi (0.7 bar)

**Outlet port:** 0.25" (6.35 mm) OD tubing – (the fitting can be easily replaced to accommodate other sizes).

### Transfer line dimensions:

Heated length about - 40" (1.0 m)  
 0.25" OD, 3/16" ID  
 Thickness w/ insulation: about 0.65"  
 (16-17 mm)

Other lengths and sizes available as options

### Transfer Line Temperature:

Maximum (Continuous) 100 °C  
 Minimum: Ambient

### Dimensions:

(W x H x D)  
 8.7" x 11.4" x 16"  
 (22.2 x 29 x 40.7 cm)

(Not including protrusions in front and back)

### Instrument Weight:

27 lb (12 kg)

### Power Requirements:

110 VAC, 400VA, 60 Hz nominal  
 (Optional): 220 VAC, 400VA, 50 Hz nominal

These specifications are subject to change at any time

Windows® is a reg. trademark of Microsoft Corporation  
 Swagelok® is a reg. trademark of Swagelok Corporation